

# Comparison of the Simplexa HSV 1 & 2 Direct Kit and Laboratory-Developed Real-time PCR Assays for Herpes Simplex Virus Detection

Jane Kuypers<sup>1,2</sup>, Gregory Boughton<sup>1</sup>, Jina Chung<sup>1</sup>, Lindsay Hussey<sup>1</sup>, Meei-Li Huang<sup>1,2</sup>, Linda Cook<sup>1,2</sup>, Keith R. Jerome<sup>1,2</sup>  
University of Washington<sup>1</sup> and Fred Hutchison Cancer Research Center<sup>2</sup>, Seattle WA, USA

## BACKGROUND

Rapid detection and differentiation of herpes simplex viruses (HSV) is important for patient management and treatment, especially in cases of HSV meningoencephalitis.

The FDA-cleared Simplexa HSV 1 & 2 Direct Assay (Focus Diagnostics, Cypress, CA) is a rapid, sample-to-result, microfluidic real-time PCR system for the direct detection and differentiation of HSV types 1 and 2 from cerebral spinal fluid specimens (CSF) that have not undergone nucleic acid extraction.

### Objective

To compare the performance of the Simplexa HSV 1 & 2 Direct assay to that of the University of Washington laboratory-developed real-time PCR tests (LDTs) for the detection of HSV1 and HSV2 DNA in CSF.

## METHODS

### CSF (n=230)

- Patients seen at University of Washington Medical Center (UW)
  - 81 residual CSF, stored at -80°C.
    - Collected March, 2009-August, 2013
    - 26 previously positive for HSV
  - 80 CSF tested prospectively
    - Collected February-April, 2014
- Neonatal herpes study (NN)
  - 69 residual CSF, stored at -80°C.
    - Collected 1981-1989
    - 34 patients
    - Median age 24 days (range, 2-258 days)
- All samples tested by Simplexa Direct and LDTs on same day.

### LDTs

- Automated nucleic acid extraction
  - Roche MagNA Pure 96 or MagNA Pure LC
  - Input vol 200 uL, elution vol 100 uL
- Real-time PCR amplification
  - Primers and TaqMan probe target HSV glycoprotein B gene
  - 7500 Real Time PCR System (Applied Biosystems)
  - 95% LOD = 250 viral copies/mL

### Simplexa Direct

- 50 uL reagent mix, 50 uL CSF
- Direct Amplification Disc
- 3M Integrated Cycle
- Primers target HSV polymerase gene
- Results ~1 hour

## RESULTS

### HSV 1 Detection by LDT and Simplexa Direct

LDT	Simplexa Direct			Total
	Pos	Neg	Invalid	
Pos	3	0	0	3
Neg	1	154	3 <sup>@</sup>	158
Total	4	154	3	161
Concordance (%)			<b>97.5</b>	
Kappa			<b>0.59</b>	
Positive agreement (%)			<b>100</b>	

<sup>@</sup>3 (1.3%) of 230 clinical CSF were invalid by Simplexa Direct.

LDT	Simplexa Direct			Total
	Pos	Neg	Invalid	
Pos	11 <sup>#</sup>	2	0	13
Neg	2	51	0	53
Total	13	53	0	66
Concordance (%)			<b>94.0</b>	
Kappa			<b>0.81</b>	
Positive agreement (%)			<b>85.0</b>	

\*3 samples pos, but not typed, by LDT. Of these, 2 were HSV2 pos and 1 was neg by Simplexa Direct.  
#8 specimens from 5 patients were positive for both HSV1 and HSV2 by LDT and Simplexa Direct.

LDT	Simplexa Direct
neg	39.2
neg	33.4
neg	35.4
38.7	neg
41.2	neg

### HSV 2 Detection by LDT and Simplexa Direct

LDT	Simplexa Direct			Total
	Pos	Neg	Invalid	
Pos	21	1	0	22
Neg	0	136	3 <sup>@</sup>	139
Total	21	137	3	161
Concordance (%)			<b>97.5</b>	
Kappa			<b>0.90</b>	
Positive agreement (%)			<b>95.5</b>	

LDT	Simplexa Direct			Total
	Pos	Neg	Invalid	
Pos	21 <sup>#</sup>	3	0	24
Neg	3	39	0	42
Total	24	42	0	66
Concordance (%)			<b>91.0</b>	
Kappa			<b>0.80</b>	
Positive agreement (%)			<b>88.0</b>	

LDT	Simplexa Direct
neg	36.4
neg	39.0
neg	40.1
33.7	neg
38.3	neg
32.7	neg
36.7	neg

### Simplexa Direct Performance

Target	Performance (%)	
	Sensitivity	Specificity
HSV1	87.5	97.2
HSV2	91.3	96.7

Target	Dil	LDT		Simplexa Direct	
		mean Ct	% pos	mean Ct	% pos
HSV1	10 <sup>-6</sup>	33.6	100	36.0	100
	10 <sup>-7</sup>	37.9	100	39.3	50
HSV2	10 <sup>-6</sup>	32.7	100	36.6	100
	10 <sup>-7</sup>	35.5	100	40.0	75

Target	N	Mean Ct	%CV
HSV1	8	22.9	1.8
HSV1	8	31.9	3.9
HSV2	8	23.6	2.4
HSV2	8	33.0	2.8
Internal control	40	30.2	0.7

## SUMMARY and CONCLUSIONS

- The Simplexa Direct performed well compared with our established LDTs for the detection of HSV1 and HSV2 in CSF.
- The majority of targets missed by both assays had high Ct values (>35), indicating low viral loads in the specimens.
- When tested using serially diluted samples, the relative sensitivity of Simplexa Direct was similar to that of LDT for HSV2. For HSV1, LDT showed slightly better sensitivity than Simplexa Direct.
- Compared to our LDTs, which require several hours for nucleic acid extraction and PCR amplification, Simplexa Direct, which provides results in ~ 1 hour, will allow for more rapid antiviral treatment and patient management.

Acknowledgement::  
Focus Diagnostics donated the Simplexa Direct HSV 1 & 2 kits for this study.